

# **FOOD MICROBIOLOGY**

### **Abstract**

This laboratory report deals with the microbial testing of Sauerkraut and Your hurt in order to identify the ***“number of mesophilic aerobes, coliforms and LAB”*** in each of the given food samples. This laboratory report also focusses on the determination of ***“percentage of lactic acid”*** produced in every given sample by using phenolphthalein as an indicator.

## **Table of Contents**

1.0 Introduction	4
2.0 Material and Method	4
3.0 Result	4
4.0 Discussion	14
References	16

## 1.0 Introduction

Sauerkraut is the most prevalent and oldest fermented foods in entire Western Europe which are mainly produced by fermentation of fresh cabbage. As fresh cabbage is the habitat of the various microbial community therefore largely contains Mesophilic aerobes and coliforms bacteria. But on fermentation, the number of these two microbial communities decreases and due to change in pH number of lactic acid bacteria (or LAB) increases in the final mixture of sauerkraut. Moreover, yoghurt also contains a large number of LAB as it is produced by the fermentation of milk product (Vasiee *et al.*, 2018). Hence, the purpose of this microbial testing is to determine the “**number of LAB, coliforms and mesophilic aerobes**” in each of the food products along with the “**amount of lactic acid**” produced in each of the samples.

## 2.0 Material and Method

### Materials needed:

NA agar, PCA agar, MacConkey agar, pH meter, phenolphthalein, spreader, Petri dish, test tube, beaker, distilled water, Bunsen burner, LAF

**Method:** The methodology for doing the entire lab report is the same as the procedure which was previously provided in the class. No changed in any of the procedure, dilution and time have been made while doing the microbiological experiments.


## 3.0 Result

### 3.1 Determination of LAB by Spread plate method

Sample	Duplicate spread Plate count at different dilution						
	$10^{-1}$	$10^{-2}$	$10^{-3}$	$10^{-4}$	$10^{-5}$	$10^{-6}$	$10^{-7}$
Sauerkraut	4	0	0				

<b>Day 0</b>	2	0	0				
<b>Sauerkraut Day 7</b>			TNTC	TNTC	327	26	1
			TNTC	TNTC	338	34	3
<b>Yoghurt</b>	0	0	65	13	1	0	0
	0	0	73	14	0	0	0

**Table 1: Number of colonies formed by LAB by Spread plate method**



**Figure 1: Colonies formed by Spread plate method in Yoghurt ( $10^{-3}$ ) dilution plate**

Therefore, for the calculation of Colony Forming Unit for LAB in each of the food samples:

$$\text{CFU} = (\text{Number of colonies in each of the plate}) / (\text{Total volume plate} \times \text{Total dilution used})$$

Food Sample	No. of Set	Dilution	Number of colonies
<b>Sauerkraut Day 0</b>	Set A	$10^{-1}$	$(40) / (1\text{ml} \times 10^{-1}) = 400 \text{ colonies}$
	Set B	$10^{-1}$	$(20) / (1\text{ml} \times 10^{-1}) = 200 \text{ colonies}$
<b>Sauerkraut Day 7</b>	Set A	$10^{-5}$	$(327) / (1\text{ml} \times 10^{-5}) = 327 \times 10^5 \text{ colonies}$
	Set B	$10^{-5}$	$(338) / (1\text{ml} \times 10^{-5}) = 338 \times 10^5 \text{ colonies}$
	Set A	$10^{-6}$	$(26) / (1\text{ml} \times 10^{-6}) = 260 \times 10^5 \text{ colonies}$
	Set B	$10^{-6}$	$(34) / (1\text{ml} \times 10^{-6}) = 340 \times 10^5 \text{ colonies}$
	Set A	$10^{-7}$	$(1) / (1\text{ml} \times 10^{-7}) = 100 \times 10^5 \text{ colonies}$
	Set B	$10^{-7}$	$(3) / (1\text{ml} \times 10^{-7}) = 300 \times 10^5 \text{ colonies}$
<b>Yoghurt</b>	Set A	$10^{-3}$	$(65) / (1\text{ml} \times 10^{-3}) = 65 \times 10^3 \text{ colonies}$
	Set B	$10^{-3}$	$(73) / (1\text{ml} \times 10^{-3}) = 73 \times 10^3 \text{ colonies}$
	Set A	$10^{-4}$	$(13) / (1\text{ml} \times 10^{-4}) =$

			130X 10 <sup>3</sup> colonies
	Set B	10 <sup>-4</sup>	(14)/(1ml X 10 <sup>-4</sup> )= 140X 10 <sup>3</sup> colonies
	Set A	10 <sup>-5</sup>	(1)/(1ml X 10 <sup>-5</sup> )= 100X 10 <sup>3</sup> colonies

**Table 2: CFU calculation for determination of the number of LAB in every sample**

### 3.2 Determination of Mesophilic Aerobe by pour plate method

Sample	Duplicate spread Plate count at different dilution						
	10 <sup>-1</sup>	10 <sup>-2</sup>	10 <sup>-3</sup>	10 <sup>-4</sup>	10 <sup>-5</sup>	10 <sup>-6</sup>	10 <sup>-7</sup>
<b>Sauerkraut Day 0</b>	142	23	0				
	129	27	0				
<b>Sauerkraut Day 7</b>			TNTC	TNTC	TNTC	329	30
			TNTC	TNTC	TNTC	343	35
<b>Yoghurt</b>	0	0	0	57	7	0	0
	0	0	0	52	3	0	0

**Table 3: Colonies formed by Mesophilic aerobes in PCA pour plate method**



**Figure 2: Colonies formed by PCA plate method in Sauerkraut Day 7 ( $10^{-7}$ ) dilution plate**

Therefore, for the calculation of colony-forming unit of mesophilic aerobes in every given food sample:

Food Sample	No. of Set	Dilution	Number of colonies
Sauerkraut Day 0	Set A	$10^{-1}$	$(142) / (1\text{ml} \times 10^{-1}) = 1420$ colonies
	Set B	$10^{-1}$	$(129) / (1\text{ml} \times 10^{-1}) = 1290$ colonies



	Set A	$10^{-2}$	$(230)/(1\text{ml} \times 10^{-1}) = 2300 \text{ colonies}$
	Set B	$10^{-2}$	$(270)/(1\text{ml} \times 10^{-1}) = 2700 \text{ colonies}$
<b>Sauerkraut Day 7</b>	Set A	$10^{-6}$	$(329)/(1\text{ml} \times 10^{-6}) = 329 \times 10^6 \text{ colonies}$
	Set B	$10^{-6}$	$(343)/(1\text{ml} \times 10^{-6}) = 343 \times 10^6 \text{ colonies}$
	Set A	$10^{-7}$	$(30)/(1\text{ml} \times 10^{-7}) = 300 \times 10^6 \text{ colonies}$
	Set B	$10^{-7}$	$(35)/(1\text{ml} \times 10^{-7}) = 350 \times 10^6 \text{ colonies}$
<b>Yoghurt</b>	Set A	$10^{-4}$	$(57)/(1\text{ml} \times 10^{-4}) = 57 \times 10^4 \text{ colonies}$
	Set B	$10^{-4}$	$(52)/(1\text{ml} \times 10^{-4}) = 52 \times 10^4 \text{ colonies}$
	Set A	$10^{-5}$	$(7)/(1\text{ml} \times 10^{-5}) = 70 \times 10^4 \text{ colonies}$
	Set B	$10^{-5}$	$(3)/(1\text{ml} \times 10^{-5}) = 30 \times 10^4 \text{ colonies}$

**Table 4: CFU calculation for determination of the number of mesophilic aerobes in every food sample**

### 3.3 Determination of Coliforms by MNP

Food Sample	Number of Coliforms
Sauerkraut Day 0	(0, 0, 0)
Sauerkraut Day 7	(0, 0, 0)
Yoghurt	(0, 0, 0)

**Table 5: Result for number of coliforms using MacConkey Broth**



**Figure 3: Result of MPN broth**

### 3.4 pH and Titration



**Figure 4: Titration by phenolphthalein**

Sample	pH value	Titre value
Sauerkraut Day 0	5.55	0.1ml
Sauerkraut Day 7	3.57	4.0ml

<b>Yoghurt</b>	4.15	6.0ml
----------------	------	-------

**Table 6: pH and % of lactic acid (titration value) for every food sample**

For the calculation of % of lactic acid in the considered food samples:

**% of Lactic acid= (“titre/ml X molarity of NaOH X mol. Mass of lactic acid”)/ (“ml of sample X10”)**

<b>Sample</b>	<b>Percentage of Lactic acid</b>
<b>Sauerkraut Day 0</b>	$(0.1 \times 0.1 \text{ M} \times 90.08) / (5 \times 10) = 0.0180 \%$
<b>Sauerkraut Day 7</b>	$(4.0 \times 0.1 \text{ M} \times 90.08) / (5 \times 10) = 0.7206 \%$
<b>Yoghurt</b>	$(6.0 \times 0.1 \text{ M} \times 90.08) / (5 \times 10) = 1.0809 \%$

**Table 7: Calculation of lactic acid in the food samples**

#### **4.0 Discussion**

The intrinsic property of Sauerkraut and Yoghurt do change the microflora in the samples. Yoghurt is a fermented product of milk and thus will contain a greater number of LAB compared to Sauerkraut which contains more aerobic bacteria. Moreover, if the fermentation was continued, the pH of Sauerkraut would have fallen more as the number of LAB would have increased in the mixture. The result of MNP clearly states that there was no coliform in the given samples. But the temperature at which the MPN method was carried out may give other bacterial growth such as Bacillus. Further, the best approach for coliform detection is the MPN method as the test provides a colour change in the MPN broth upon coliform detection. Hence, MPN provides advantages over the traditional plate method as by MPN coliform bacteria can be easily detected and differentiated from other aerobic bacteria growing at 37°C. The pH and titre value

ratio of both the sample was almost similar as both the samples contained LAB as the primary microbial community.

## References

- Alpaslan, E., Geilich, B.M., Yazici, H. and Webster, T.J., 2017. pH-controlled cerium oxide nanoparticle inhibition of both gram-positive and gram-negative bacteria growth. *Scientific reports*, 7(1), pp.1-12.
- Khan, M.I., Mukherjee, K., Shoukat, R. and Dong, H., 2017. A review on pH sensitive materials for sensors and detection methods. *Microsystem Technologies*, 23(10), pp.4391-4404.
- Lee, S., Lee, S., Singh, D., Oh, J.Y., Jeon, E.J., Ryu, H.S., Lee, D.W., Kim, B.S. and Lee, C.H., 2017. Comparative evaluation of microbial diversity and metabolite profiles in doenjang, a fermented soybean paste, during the two different industrial manufacturing processes. *Food chemistry*, 221, pp.1578-1586.
- Taylor, B.C., Lejzerowicz, F., Poirel, M., Shaffer, J.P., Jiang, L., Aksenov, A., Litwin, N., Humphrey, G., Martino, C., Miller-Montgomery, S. and Dorrestein, P.C., 2020. Consumption of fermented foods is associated with systematic differences in the gut microbiome and metabolome. *Msystems*, 5(2).
- Vasiee, A., Behbahani, B.A., Yazdi, F.T., Mortazavi, S.A. and Noorbakhsh, H., 2018. Diversity and probiotic potential of lactic acid bacteria isolated from horreh, a traditional Iranian fermented food. *Probiotics and antimicrobial proteins*, 10(2), pp.258-268.
- Wuyts, S., Van Beeck, W., Oerlemans, E.F., Wittouck, S., Claes, I.J., De Boeck, I., Weckx, S., Lievens, B., De Vuyst, L. and Lebeer, S., 2018. Carrot juice fermentations as man-made microbial ecosystems dominated by lactic acid bacteria. *Applied and environmental microbiology*, 84(12).